



CRISTAL Pro[®]

THE NEW
360°
BODY
CONTOURING
METHOD

A UNIQUE CONCEPT



4 simultaneous applicators



Integrated ultrasound



Body & face platform

by **Deleo**
CREATOR IN
AESTHETIC MEDICINE

CE
2797

OFFER YOUR PATIENTS THE EXPERTISE OF 360° BODY CONTOURING



All body & face areas



For women and men



At any age



All year round



Dr. SCEMAMA
Cosmetic doctor,
Paris, FRANCE

« Trusting Deleo's laboratory since the beginning, and fully satisfied with the 1st CRISTAL, I was seduced by the new technology and innovations of the CRISTAL Pro : the number of handpieces allowing to treat any fat fold, the performance and safety of its handpieces, the possibility to treat up to 4 areas at the same time, ... »



REVOLUTIONISE YOUR PRACTICE

CRISTAL Pro IS...



Dr. SARFATI
Ophthalmologist,
Toulon, FRANCE

« This machine has broken new ground in the treatment of adipocyte apoptosis. Indeed, it offers global patient management : by an ultrasound analyse of the areas to be treated; by the creation of numerous handpieces adapted to more or less accessible body areas ; by more precise management of lowering temperature, which becomes safe and scalable. »



4 SIMULTANEOUS APPLICATORS

- Time saving: treatment of up to 4 body & face areas
- New «Soft Comfort» designed applicators for fewer side effects



OPTIMIZED EFFECTIVENESS

- Homogeneous cooling down to -12°C thanks to patented CCDSTM technology
- Intelligent suction adjustment



SOFTWARE INTELLIGENCE

- Customized configuration
- Connected technology
- Ergonomic & intuitive software



360° BODY & FACE APPLICATORS

CRISTAL Pro offers a wide range of morphological applicators adapted to large areas (abdomen, love handles, thighs,...) but also to small areas (chin, arms,...).

The goal? To offer 360° coverage of the silhouette for optimized results and maximum patient satisfaction.

AGATE



Very small fat folds for face

Indication : Chin

JADE



Small to medium fat folds for body

Indications : Bra folds, arms, knees, calves

SAPPHIRE



Medium to large fat folds

Indications : Abdomen, love handles, inner thighs

AMETHYST



Applicator without suction, specifically dedicated for flat areas

Indications : Arms, saddle bags, inner and outer thighs, calves

AMBER



Small fat folds for body

Indications : Lower jaw, pectoral area, knees

EMERALD



Medium fat folds

Indications : Bra folds, arms, small love handles, knees

RUBY



Large fat folds

Indication : Abdomen



Dr. ADDA

Cosmetic doctor,
Boulogne-Billancourt, FRANCE

« As a CRISTAL Pro user, I am delighted with the technical features of this machine, both ergonomically and functionally. The first results are convincing. »




EXCLUSIVE INTEGRATED ULTRASOUND SCANNER

The CRISTAL Pro ultrasound scanner allows for the precise measurement of adipose tissue before and after a treatment.

This linear dynamic probe, emitting at a central frequency of 7.5 MHz, can be used for the visualization of different skin layers, including adipose tissue, but also for measuring muscle volume.

The CRISTAL Pro ultrasound scanner is simple and easy to use.





A UNIQUE CONCEPT

CRISTAL Pro accompanies you at every stage of your practice development.
Joining the CRISTAL Pro community guarantees personalized and long-lasting support.



TRAINING

In-depth initial training for your entire team

Continuing training thanks to our committee of scientific experts



AFTER-SALES SERVICE

Dedicated hotline

Quick intervention & spare parts' shipment within 48 hours



PREMIUM MARKETING PACK

All the necessary elements to communicate in your clinic

Media and digital coverage (TV, press, display, social networks, blog)



PATIENT RECRUITMENT SUPPORT

Professional space on www.partnerzone-deleo.fr to download all technical, clinical & marketing support material

BODY CONTOURING, YOUR SUCCESS VECTOR




By 2022, the Body Contouring market will grow by 14.7% per year¹



Cryolipolysis stands out with an unbeatable growth of 26.6%¹

2x
80.6% of CRISTAL patients wish to have a second session³

9 million Body Contouring treatments performed this year²



“

Dr. ESSAYAGH

Cosmetic doctor,
Le Cap d'Antibes, FRANCE

« We are pleased to be equipped with the Cristal Pro cryolipolysis system to deliver on the promise of reducing localized fat deposits in a serious, effective and professional way. The ultrasound module allows us to refine the diagnosis of the right indications and thanks to the precise measurement of suction pressures that allows us to choose the right handpiece, the right cooling temperature, all in complete safety. »

”

1 Source: Aesthetic Insight – August 28, 2018
2 Source: Data Tribune 2017 (MII, DRG1) Non-invasive Fat reduction, Research & Market, Aug 16
3 Source: Michaël Naouri M. Fat Removal Using a French Cryolipolysis Device: Retrospective Study of 418 Treatments. France. 2016.

CRISTAL Pro, THE STARTING
POINT OF THE **CRISTAL**
BODY HEALTH CARE PROGRAM

1 REDEFINE YOUR SILHOUETTE
with CRISTAL Pro®, 360° Body Contouring

Crystallizes and destroys adipocytes.
www.cryolipolysecristal.com

2 SCULPT YOUR SILHOUETTE
with CRISTAL Fit®, the High Intensity
Electromagnetic Stimulation Technology

Tones and shapes muscles.
www.cristal-fit.fr/en

3 TIGHTEN YOUR SKIN
with DIAMOND®, combining infrared thermal
energy and modular suction

Tightens the skin and reduces cellulite.
www.diamond-cellulite.com

4 PERPETUATE AND POTENTIATE
YOUR RESULTS

with a healthy lifestyle combining :

- Regular physical activity;
- A balanced diet;
- An adapted care routine.

Discover advice from our experts in the CRISTAL®
post-treatment guide.



VISIBLE RESULTS WITH CRISTAL TECHNOLOGY

Choosing a technology from the CRISTAL range insures visible results after just a few sessions, maximizing your patients' satisfaction.

@DR.SCHINAZI



Before

ABDOMEN

After 1 session

@DR.SCEMAMA



Before

**LOVE
HANDLES**

After 1 session

@DR.MALET




Before

SADDLEBAGS

After 1 session

@DR.SARFATI



Before

**DOUBLE
CHIN**

After 1 session



TECHNICAL SPECIFICATIONS



Dr. AZOULAY

Cosmetic doctor,
Boulogne-Billancourt, FRANCE

« After 20 years of expertise in aesthetic medicine, it seemed logical to me to equip myself with the best in terms of cryolipolysis technology, so I chose the CRISTAL Pro, which is the most advanced. »



Temperature

0°C to -12°C
(1°C interval)

Number of applicators available

7

Number of simultaneous applicators

4 with
bi-zone
control

Cooling circuit

Mono
Ethylene
Glycol

Power

1200 W

Pressure

0 mbar to
-700 mbar
(10 mbar
interval)

Power supply

AC 110-240V
~50/60Hz
16-8A

Software

Large 15.6-inch
screen
Android PC
7.1.2 tactile



CRISTAL Pro[®]

With more than 600 installations in France and around the world and more than 1 000 000 treatments performed, CRISTAL is the reference and the European leader in Body Contouring.

To bring a global treatment solution to patients, CRISTAL evolves in its Pro version and accompanies you towards your practice's success.

www.deleo-medical.com



Dr. DELLIERE

Plastic surgeon,
Cesson-Sévigné, FRANCE

« We are very satisfied with our CRISTAL cryolipolysis device in our centre.

The CRISTAL Pro allows us to go even further in terms of treatment zones with its new applicators, in comfort for the patients and in efficiency with the possibility to use 4 applicators on the same session! »



Deleo
CREATOR IN
AESTHETIC MEDICINE

Deleo S.A.S.
239 via Nova
Pôle d'Excellence Jean Louis
83600 Fréjus, FRANCE
T. 04 94 45 83 75 | contact@deleo.fr
www.deleo-medical.com

